GLOBUS

GENESY 3000

GRAPHIC BACKLIGHT DISPLAY:

Genesy 3000: 67

Rehabilitation, 35 Pain, 34 MENS, 9 Denervated muscles, 29 Kotz, 25 Interferential, 1 Ionophoresis, 100 Training, 64 Fitness, 119 Beauty

DENERVATED DIAGNOSTIC

FUNCTIONAL RE-EDUCATION: the

programs are aimed at correctly carry out the rehabilitation of the injured person by means of specific programs having rectangular, kotz, denervated and interferential waves.

IONOPHORESIS (Galvanic Current) : IONOPHORESIS is suitable for the treatment of inflammatory states like bursitis, tendinitis, fascitis, myositis etc.

MICRO-CURRENT (MENS): Genesy 3000 presents the new and innovative micro-current, with their regenerative and curative effects, particularly indicated for Sports, Beauty and above all Health applications.

PROGRAMMABLE: 15 empty memories.

POWER SUPPLY BOTH BY RECHARGEABLE BATTERIES AND CONNECTION TO THE CURRENT NETWORK.

2+2 FUNCTION: possibility of carrying out 2 different programs simultaneously.

MULTI USER: Genesy 3000 allows the sharing of the electrostimulator among 30 users, keeping the respective setting, programs and preferences unchanged.

AUTOMATIC WAY (AUTO STIM®): Genesy 3000 is the only electrostimulator that allows you to carry out a program in automatic way, without any manual intervention on the intensity regulation.

MEMO FUNCTION: Genesy 3000 allows you to create a library of personalized programs.

SYNCROSTIM: Possibility of increasing the intensity of the 4 channels simultaneously by pushing an unique pushbutton.

LAST 10: Genesy 3000 keeps in memory the last executed programs. In this way they can be carried out again in a simple and fast way.

FULL RANGE OF PRODUCTS: Genesy 3000 has a series of programs addressed to the functional recovery (atrophy, force and mobility) and to the general well being (TENS and MENS) and lonophoresis. It also includes treatments for sport and beauty.

PROGRAMS SPORT AND FITNESS: These programs complete the main range of programs for the development of the force (maximum force, explosive force, resistance etc.) and specific beauty treatments (lipolyse, drainage, tonification etc.).

GLOBUS

PROGRAM LIST

REHABILITATION PROGRAMS

Atrophy recovery Motor point pen Muscle mass recovery Strengthening Mobilization Agonist/Antagonist Cruciate ligament elesticity Achille's tendon elasticity Upper limbs – ligaments elasticity Upper limbs - hemiplegia Lower limbs - hemiplegia Vescical instability Mixed incontinency Stress incontinency Paraplegia (wound prevention) Lower limbs spasticity Upper limbs spasticity Hip prosthesis

PAIN PROGRAMS

Anti-accomodation analgesic tens Fracture pain Algodistrophy Fantom limb 1 and 2 Cronic pain Rib fracture Trapezius pain Hip pain Knee pain Ankle pain Endorphinic Tens Haematoma Contracture Sciatica Antistress Lumbago Cervico-brachial neuralgia Cervicals Scapulohumeral periarthritis Epicondylitis Nausea Carpal tunnel Hip arthrosis Rotula chondritis Menstrual pain Intense menstrual pain Acute pain analgesic tens Chronic pain analgesic tens Modulated Tens Backbone compression Periostitis Finger pain Rizopathy

IONOPHORESIS lonophoresis

MICRO-CURRENT PROGRAMS Muscular contracture

Muscular strain Muscular sprain Muscular energetic restoration Contusion Rheumatics Arthritis Oedema Herpes Zoster Hematoma Ischemic ulcer of the skin Sciatica Femoral neuralgia Lumbago Slipped disc and lumbar spondylosis Brachial neuralgia Stiff neck Migraine "Whiplash" neckpain Cervical spondylosis Shoulder sprain Scapulohumeral periarthritis Epicondylitis Carpal tunnel syndrome Hand arthritis Finger fracture Brest micro-lifting Decolleté micro-lifting Face micro-lifting Face skin elasticity Knee sprain Rotula tendon inflammation Knee osteoarthritis Ankle sprain Achille's tendon inflammation

Demo

Rectangular pulse waveform pre set program Rectangular 1 - 2 - 3Triangular pulse waveform pre set program Triangular 1-2-3Trapezoidal pulse waveform pre set program Trapezoidal 1-2-3

KOTZ Muscular training Coxarthrosis Fatigue recovery Gluteus deficit COLLES fracture and flexors reactivation Toracic hyperkyphosis Muscular hypertrophy Dorsal para vertebral muscular hypotrophy L5 pain residual hypotrophy Quadriceps anterior rectum hypotrophy Post-trauma ischium-crural hypotrophy Deltoid post-trauma hypotrophy Triceps hypotrophy Vastus medialis hypotrophy Lumbago Not surgery medial meniscus pain Vertebral osteoporosis Flat foot Lumbar-thoracic fracture after-effects Quadriceps Muscular strengthening after meniscus pain Metatarsus pain Gluteus reactivation Serratas muscular strengthening Quadriceps isometric strengthening Warm up Scoliosis Muscular toning

INTERFERENTIAL PROGRAMS

Steady (for painful tendon and muscular pathologies) Sweep 80 - 100 Hz (for painful tendon and muscular pathologies) acute phase Sweep 5 - 100 Hz (for painful tendon and muscular pathologies) acute phase and inflammatory state Sweep 5 - 20 Hz (for painful states where edemas reabsorption is required) nump effect Multi freq (edemas reabsorption in the presence of acute pain) Sweep 5 - 50 Hz (for muscular pathologies where edemas reabsorption is needed) Sweep 50 - 100 Hz (chronic and acute muscular pain) LP (used as alternative to the diadynamic LP) CP (used as alternative to the dynamic COST, LP, CP (Association of the precedine programs) AMF (recommended to treat wide painful areas) Knee degenerative arthrosis Sub-acromial bursitis Myotensive headhache Acute cervical-brachial pain Coxarthrosis Knee sprain Acute cervical pain Acute dorsal pain Acute epycondilitis Chronic lumbago Calcific scapulohumeral periarthritis Contractures Shoulder sprain Hankle sprain

SPORT PROGRAMS

Demo Capillarisation Warm up Pre-competition warm up Active recovery Maximum force Resistance force Explosive force Aerobic resistance Reactivity Post-competition recovery Relaxation Ipertrophia Dynamic work

FITNESS PROGRAMS

Firming Sculpting Toning Mass increasing Body Sculpting Body definition Jogging Anaerobic fitness Aerobic fitness Cramp prevention

BEAUTY PROGRAMS Drainage

Lipolysis Relaxant massage Toning massage Energetic massage Connective massage Vascolarisation Post-pregnancy drainage Post-pregnancy lipolysis Post-pregnancy firming Breast firming Breast firming Swallen arms Face capillarisation Face lifting-effect Arterial failure Venous failure Definition Body shape

TECHNICAL FEATURES

Channels: 4 indipendent channels (8 electrodes) Frequency: 0,3 - 150 Hz Power: 120 mA per channel Amplitude: $45 - 450 \ \mu s$ Wave form: symmetrical biphasic square pulse Power supply: battery charger from 18 V -120 VA

Battery: 12 V. 2000 mAh Conformity : the device complies with the ISO 9001, EN 46001 standards for electromedical products that are MDD (Medical Device Directive) certified.

ACCESSORIES :

- 1 bag 1 GENESY 3000 electrostimulator 4 electrode connection cables 1 battery charter 4 reusable adesive electrodes 50 x 50
- 4 reusable adesive electrodes 50 x 90
- 1 operating manual

CURRENT TYPES

- EMS TENS
- MCR
- KOTZ
- DENERVATED
- INTERPHERENTIAL IONOPHORESIS

DENERVATED PROGRAMS