ENDOWMENT

1x — Unit BE-28T

- 2x Wires for electrodes connection
- 4x Self-adhering electrodes, with
- 2 mm socket connection 2x — Alkaline batteries AA1.5 V (LR6)
- 1x Instructions booklet
- 1x Stuffed bag with compartments



WARNINGS!

BE-28T is a medical device CE. Carefully read the warnings and the instructions for use. The unit cannot be used: - on patients with cardiac pacemakers; - on pregnant women - on areas affected by inflammatory processes.

TECHNICAL SPECIFICATIONS

Supply:2 alkaline batteries1,5 V AA (LR6)

Output:Maximum 130mApp across 500 Ohm load

Pulse Width: From 50µS to 250µS adjustable in steps of 10µS

Pulse rate: From 1Hz to 150Hz (1-2-3-4-5-10-12-14-16-18-20-25-30-35-40-45-50-60-70-80-90-100-110-120-130-140-150)

Waveform: Symmetrical Bi-phasic rectangular - Asymmetrical Biphasic rectangular - Mono-phasic rectangular

Treatment Timer: Continuous, 15min, 30min, 45min, 60min, 90min.

Dimensions: 138mm x 6 8mm x 28mm

Weigth: 160 g. including batteries



IntelliSTIN BE-28T

0051

Professional unit for pain management through electrical stimulation (T.E.N.S.).

Since 1981...Hardware & Software for Rehabilitation

We exactly know what you expect by a professional analgesic unit...

...and we developed...

By a modern analgesic electrical stimulation unit (T.E.N.S.) a professional user expects that it is: **effective** (for achieving therapeutic aims), **easy to operate** (for avoiding the user to waste his time), **flexible** (for allowing to treat any pain situation), **robust** (to resist for long time). Further to the above, **IntelliSTIM® BE-28T** also offers:

- 6 Large LCD display
- 6 Intuitive Operation
- 6 9 operating programs (editable)
- 6 any parameters shown
- 6 "edit & feel" function
- 6 ACUTE and CHRONIC pain protocols
- 6 3 Selectable wave shapes





Professional unit for pain management.

- Dual intensity-independent channels
- Last used program loading
- Count-down timer with switch-off at the end
- Auto-switch-off if not in use
- Accumulative recording of the treatment time (in hours)
- Electrodes protection circuit (just in the model BE-28TC)
- Lock function

ENDORPHINS or "Gate –Control"?

Among the other parameters, on the display, it is also shown the applied analgesic technique, according to the setup of the unit.