#### **AVAILABLE MODELS**

#### ■ MAC 1254 - MAGNETOMED 2000

Generator for magnetic therapy assembled on table container with two independent outputs.

# **MAGNETOMED**

#### **TECHNICAL FEATURES**

- 20 stored programs
- 10 storable programs on the internal memory
- Maximum power per each output: 100 Gauss
- Adjustable frequency from 1 to 100 Hz
- Adjustable treatment time: 99 minutes

## **MAGNETIC THERAPY**

## STANDARD ACCESSORIES

Power cable	I
User manual	l
Fuses	2 da 315mA-T
Ring Magnet	l
(for the device's test)	

## ACCESSORIES TO BE CONNECTED TO THE GENERATOR

STATIVE APPLICATOR EQUIPPED WITH VANE SLIDING SOLENOID, PERFECT FOR TREATMENTS WHILE BEING

• ACC404 - SA





# STATIVE APPLICATOR EQUIPPED WITH 50 CM CYLINDER

• ACC598 - SC50



# 3 SHELVES TROLLEY

• ACC604




#### **MAGNETOMED - OPTIONAL ACCESSORIES**

High-Tech aluminium bed with 2 highly sliding rails and 2 cylinders (a cylinder with a 60 cm diameter and 36 cm width and the other with a 60 cm diameter and 22 cm width) imitation leather lined and provided with the exclusive MFC patent concentrating the lines of force inside the solenoid without wasting them outside (in compliance with the legislative decree 626). The North Pole is to be found for a correct orientation...

• ACC544 - LC60/2

High-Tech aluminium bed with I highly sliding rail and I cylinder with a 60 cm diameter and 36 cm width, imitation leather lined, provided with the exclusive MFC patent concentrating the lines of force inside the solenoid without wasting them outside (in compliance with the legislative decree 626). The 60 cm cylinder lies on ball bearings and can slide 100% of the total run. The North Pole is to be found for a correct orientation.





• ACC543 - LC60

Portable cylinder, 30 cm diameter, equipped with foot for level surfaces and provided with exclusive MFC patent concentrating the lines of force inside the solenoid without wasting them outside (in compliance with the legislative decree 626). The North Pole is to be found for a correct orientation.

• ACC549 - CP30



Portable cylinder with 50 cm diameter for lower and upper limbs treatment. Portable cylinder, 50 cm diameter, equipped with foot for level surfaces and provided with exclusive MFC patent concentrating the lines of force inside the solenoid without wasting them outside (in compliance with the legislative decree 626). The North Pole is to be found for a correct orientation.

• ACC550 - CP50





Two applicators,  $16 \times 10 \times 3.5$  for localized treatments. The two applicators can either work together (for big areas) or in opposition (to focus the magnetic induction into a specific area). Plastic container, lined with grey imitation leather to find the North Pole and blue for the South Pole. Light weight and soft colour while in action.

• ACC615 - CP



