

VACUUMED

DEDICATED TO
PEOPLE LOOKING
FOR THE TOP

The Vacuum therapy is a therapeutical method which uses the effect of the vacuum in the body by combining the improvement of the tissue microcirculation with an anti-inflammatory and detoxing action.

Concretely the vacuum therapy consists in applying an inner tube to a body part, in which a depression is created and afterwards a compression to stimulate the venous and lymphatic circulation: this causes a suction effect which increases the blood afflux in the involved area, stimulating the circulation and the oxygenation, and accelerating the turnover of substances and the elimination of toxins.

In the peripheral tissues the oxygen levels increase substantially and the carbon dioxide decreases: by improving the blood flow the exchange activity in the capillaries is favoured and by accelerating the blood flow the venous stasis which is the main cause of varicose veins and venous insufficiency.

The vacuum application combined with the stimulation currents will be carried out using vacuum electrodes. Such electrodes are connected through pipes to a suction pump which causes a negative pressure inside the electrodes. During this action, the high depression requires a great quantity of blood (hyperaemia) thereby improving the therapy effect through the different stimulation currents.

INDICATIONS

Such therapy is clearly useful for treating all forms of vascular diseases affecting the legs: arteriosclerosis of lower limbs, Burger, Raynaud; and all forms of venous and/or lymphatic insufficiency, such as swollen legs, lymphedema, phleboedema, and cellulitis.

Vacuumed can be therefore connected to the devices for the electrotherapy:

- THERAPIC 7200 - 7400 / 8200 - 8400 / 9200 - 9400
- THERAPIC 1500 - 2000
- ELEKTRA 2 - 4
- COMBIMED 200 - 2200



• MAC516 - VACUUMED

DEDICATED TO
PEOPLE LOOKING
FOR THE TOP



TECHNICAL FEATURES

- Suction from 0 - 0,6 bar
- Frequency from 10 a 100 pulses/minutes
- Power supply 230V 50/60Hz
- Absorption 40VA

SUPPLIED ACCESSORIES

VACUUMED

Mains cable	1
User manual	1
Vacuum cups Ø 60 mm	4
Cup sponges Ø 60 mm	8
Output cables	4
Electrotherapy connection cable	1

VACUMED - OPTIONAL ACCESSORIES

VECUUM CUPS

- ACC913 - 30 mm
- ACC695 - 60 mm
- ACC914 - 90 mm



VACUUM SPONGES

- ACC915 - 30 mm
- ACC621 - 60 mm
- ACC916 - 90 mm



VACUUM OUTPUT CABLE

- ACC696



SHORT VACUUM CABLE 20 cm

- ACC699



CONNECTION CABLE WITH ELECTROTHERAPY

- ACC697



VOCAL KIT

- ACC501

