

# EasySTIM®

## NMS -28




**2 channels low-cost Neuro-Muscular electrical stimulation unit, for professional or personal use.**

- o *Dual independent channels*
- o *Last used program loading*
- o *Count-down timer with switch-off at the end*
- o *Auto-switch-off if not in use*
- o *Accumulative recording of the treatment time (in hours)*
- o *Electrodes protection circuit (model NMS-28C only)*
- o *Lock function capability, by the doctor.*

### Technical Features

<b>Supply</b>	2 alkaline batteries 1,5 V AA (LR6)
<b>Output</b>	Max 100 mA <sub>app</sub> over 1000 Ohm load (using 200 µS pulse width)
<b>Pulse width</b>	From 50µS to 400µS adjustable in steps of 50 µS
<b>Frequency:</b>	From 1 Hz to 110Hz (1-2-3-4-5-10-12-14-16-18-20-25-30-35-40-45-50-60-70-80-90-100-110)
<b>Wave shape:</b>	Biphasic Simmetrical Pulses Monophasic Alternated Pulses
<b>Treatment timer:</b>	Continuous, 15min, 30min, 45min, 60min, 90min selectable
<b>ACTION/ REST time:</b>	From 1 to 40 sec, adjustable in steps of 1 sec.
<b>RAMP time:</b>	From 1 to 5 sec, adjustable in steps of 1 sec.
<b>Dimensions:</b>	138mm x 68mm x 28mm
<b>Weight:</b>	160 g. including batteries



### Endowment

**EasySTIM NMS-28 is supplied complete with:**

- 1x Therapeutical unit NMS-28**
- 2x Couple of wires for connecting the electrodes**
- 4x Self-adhering electrodes 50x50mm with 2mm socket connections**
- 2x Alkaline batteries AA1.5Volt (LR6)**
- 1x Instructions booklet**



Buy



Business



Related



**BEAC Biomedical**

Via Monte Bianco 12 - 27040 Portalbera (PAVIA) Italy

EN ISO 9001--EN 46001 Certified Quality Management System



For obtaining further information: Ph.: +39-(0)385-43378 Fax: +39-(0)385-48923 E-mail: info@beac.it

