

ACTIVA 500

• **MY TRAINER:** shows the images of the electrodes' placement directly on the screen.


• **389 PROGRAMS:** 13 programs G-Pulse ©, 6 Face, 37 Vascular, 228 Fitness, 95 Beauty/Aesthetics, 10 Tens.

• MICROCURRENTS

(AESTHETICS): Activa 500 presents the new and innovative "regenerating" G-Pulse © currents proposed for Beauty / aesthetics applications.

• **FACE AREA:** 6 programs developed by the CRG to help tone your face.

• **ELECTRODE HEAD:** Included a special electrode head specially designed for the micro-currents face applications.



• GRAPHIC LCD DISPLAY

• **FUNCTION 2+2:** possibility of carrying out 2 different programs simultaneously

• **MULTI-USER:** Activa 500 allows to share the electrostimulator among different users, keeping the respective setting, programs and preferences saved on the unit.

• **AUTOMATIC MODE (AUTO STIM®):** Globus Activa 500 is the only stimulator that allows to carry out a program automatically without the need of a manual intervention on the intensity regulation.

• **FUNCTION MEMO:** Activa 500 allows you to create a library with your personalized programs.

• **SYNCROSTIM:** with this function, simply by pressing one button, the intensity of all the channels increases.

• **LAST 10:** Activa 500 keeps in memory the last executed programs. In this way they can be carried out again in a very simple and fast way.

• RECHARGEABLE BATTERIES

home fitness

PROGRAMS LIST :**PHYSICAL FITNESS PROGRAMS
TONIFICATION**

Demo
Motor point pen
Active recovery
Maximum force
Relaxation
Firming 1-2
Firming Bio Pulse 1-2
Sculpting 1- 2
Sculpting Bio Pulse 1- 2
Tonification 1-2
Body Building 1-2
Body Sculpting 1-2
Definition 1-2
Jogging
Muscular initiation
Anaerobic fitness
Aerobic fitness
Cramp prevention
Aerobic resistance

VASCULAR PROGRAMS

Drainage 1-2
Drainage Bio Pulse
Drainage post- pregnancy
Swollen arms
Arterial failure
Venous failure
Vascularisation 1-2
Capillarisation

BEAUTY

Lipolysis 1-2
Lipolysis post- pregnancy
Body shape 1-2
Relaxant massage
Bio Pulse massage
Toning massage 1-2
Energetic massage 1-2
Connective massage 1-2
Post-pregnancy firming
Definition 1-2

FACE / BREAST PROGRAMS

Firming Breast
Modeling breast
Capillarisation face
LIFTING face

PROGRAMMI DOLORE

Endorphinic Tens
Acute pain analgesic Tens
Chronic pain analgesic Tens
Modulated Tens

Cervicals
Menstrual pain
Hip pain
Knee pain
Ankle pain
Cronic pain

G-PULSE PROGRAMS

Micro-lifting breast
Micro-lifting décolleté
Micro-lifting face
Skin elasticity
Cellulite
Collagene-Bio Skin
Wrinkles
Eyes wrinkles
Expression wrinkles
Stretch marks
Nails strengthening
Bio Peeling
Haematoma

TECHNICAL FEATURES

Channels: 4 independent channels (8 electrodes)
Frequency: 0,3 – 150 Hz
Power: 120 mA/channel
Impulse amplitude: 45 – 450 µs
Wave form: Symmetrical biphasic square pulsed
Power Suppli: Battery charger

Conformity: The device complies with the ISO 9001, EN 46001 standards for electro medical devices that are MDD certified (Medical Device Directive).

SUPPLY WITH :

1 transport bag
1 Electrostimulator unit ACTIVA 500
4 Electrode Cables
2 Grey Cables for microcurrents
1 Electrode Head for microcurrent
4 square electrodes
4 facial electrodes
4 rectangular Electrodes
1 Gel
1 Battery charger
1 User manual