

GLOBUS

EASY TENS

ELIMINATING THE PAIN – TENS HEALTH WELL-BEING



ELIMINATING THE PAIN

The electro stimulation is of primary importance in the elimination of the pain of articular, muscular or nervous origin. The TENS currents in fact stop quickly the pain transmission, promoting the endorphins release, that carry out a strong antalgic action that lasts in the time.



BACKHACKE AND SCIATICA

Immediate relief of the lumbar pain joined to a beneficial relaxing action of the sore muscles.

PERYARTHROSIS: SHOULDER, ELBOW AND KNEE

Elimination of the pain and of articular rigidity of this painful and inflammatory pathologies.



NECK PAIN

Relaxant action localized on the cervical area of the spinal column that eliminates the muscular spasm, the neck pains, the neck rigidity and the headaches.

LYMPHATIC DRAINAGE

The use of draining currents removes the intra-cellular liquids in excess and helps solving the circulation disorders by eliminating the limbs swelling and heaviness.



POST-TRAUMA MUSCULAR RECOVERY

All muscular system traumas and articular pathologies always involve an important loss of the mass and of the tone of the injured muscles. Thanks to EMS currents you could quickly rehabilitate the muscular function you have lost, recovering the movement ability completely.

A PRECIOUS ALLY FOR YOUR HEALTH :EASY TENS is the ideal electro stimulator for protecting your health and giving back well-being to your body, in a comfortable and safe way.

IMMEDIATE RELIEF: TENS currents eliminate the pain quickly and, day by day they help you to recover energy and vitality.

EASY AND POWERFUL: The big backlight display and the handy and intuitive controls make it very easy to use since the first treatment you carry out.

ALWAYS WITH YOU: EASY TENS is so light and compact that will be the indispensable partner of your vitality and well-being in every moment of the day.

THE BEST PROGRAMS AGAINST THE PAIN

- 1 TENS
- 2 TENS ENDORPHINIC
- 3 SCIATICA
- 4 BACKACHE
- 5 NECK PAIN
- 6 PERYARTHRTIS
- 7 DRAINAGE
- 8 MUSCULAR RECOVERY



ACCESSORIES

- PALM ELECTROSTIMULATOR
- BAG FOR THE TRANSPORT
- 2 CONNECTION CABLES
- 4 REUSABLE DHESIVE ELECTRODES 45X45 MM
- 4 REUSABLE DHESIVE ELECTRODES 45X90 MM
- 3 1,5 V AAA BATTERIES
- OPERATING MANUAL

TECHNICAL FEATURES

Channels	2
Intensity	100 lpp
Wave form	Symmetrical biphasic square pulse
Frequency	1— 90 Hz
Amplitude	50 –200 μ s
Power supply	3 1,5 V AAA batteries type LRO3
Conformity	The device complies with the ISO 9001, EN 46001 standards for electro medical devices that are MDD (Medical Device Directive) certified.

 **medizin**
comércio de equipamentos
médicos e hospitalares, lda
Rua Artur Napoleão, 226
4460-246 Senhora da Hora
www.medizin.pt • geral@medizin.pt
TEI 229 563 100 • FAX 229 563 102